



**PROVEN ON THE BATTLEFIELD.
TRUSTED IN THE CATH LAB.**



**Visit us at ACC at
booth #1659 and at
SIR at booth #2307
for a free sample.**

**For more information,
call 203-294-0000 or visit
www.z-medica.com.**

No matter the uniform — or the conditions, doctors everywhere strive to give their patients the best possible care.

That's why medical professionals all over the world trust QuikClot® Interventional™ Hemostatic Bandage to stop the bleeding fast. Proven in even the toughest conditions, the active ingredient in QuikClot® is kaolin — an all-natural, non-allergenic mineral — that activates factors in the blood to initiate and promote rapid coagulation in as little as five minutes.*

QuikClot® is the hemostatic agent of choice used by the US Military, and has been cleared by FDA for external use for moderate to traumatic bleeding. Trust QuikClot® in your Cath Lab.

- ★ **100% success in controlling hemorrhage following percutaneous vascular access.***
- ★ **Achieves hemostasis in 5 minutes or less in both interventional and diagnostic procedures using up to 12F.***
- ★ **Safe to use and non-allergenic — contains no animal, human, or botanical proteins.**
- ★ **Triggers the intrinsic pathway of the body's natural clotting cascade.**

Z-MEDICA

Z-Medica Corporation | 4 Fairfield Boulevard | Wallingford, CT 06492

QuikClot®

THE BLEEDING STOPS HERE.